


**Uncompromised
Hemodynamics.
By Design.**

A collaboration



Overview

Hemodynamics



Small Annuli



**Coronary
Access**



Durability



**Bespoke
Delivery System**





Uncompromised Hemodynamics. By Design.

ALLEGRA™ has been designed for uncompromised hemodynamic performance with single digit mean pressure gradients and high effective orifice areas¹⁻³.



unique
self-expanding
stent frame




The ALLEGRA™ has a unique self-expanding stent frame with convex and concave areas that provide high EOA (effective orifice areas), even in small annuli¹⁻³. This unique stent shape also means that there is no interaction with the ascending aorta.




The ALLEGRA™ stent design incorporates tailored radial force distribution with high radial force in the sealing area to achieve secure anchoring and low radial force in the outflow section.

compression by native annulus



wider valve opening



The low radial force in the out-flow section supports wider valve opening and contributes to the large EOA and low mean pressure gradient¹⁻³.



The ALLEGRA™ is a supra-annular valve. The new valve plane sits above the constricted and diseased native valve. This supra-annular valve position also contributes to the large EOA and low mean pressure gradient¹⁻³.



supra-annular
valve position



Maximizing the Effect in Small Annuli.



The outstanding hemodynamics of the **ALLEGRA™** are particularly important in small native annuli and surgical valves (ViV).



Patients with small annuli derive particular benefit from the use of the ALLEGRA™ valve, due to the combination of wide valve opening and supra-annular valve plane.

The ALLEGRA™ can be used with annular diameters as low as 19 mm in native valves and 16.5 mm in surgical valves. Even in these very small annuli, low single digit mean pressure gradients and high effective orifice areas can be achieved^{1,3}.



ALLEGRA™ shows better in-vitro hemodynamic performance than other TAVI devices in the Valve-in-Valve setting.


In-vitro comparisons have shown higher EOA compared to other intra-annular and supra-annular valves^{4,5}.



┌

The design features of the ALLEGRA™ make it an excellent choice for Valve-in-Valve procedures with low frame height and unique visualization of the new valve plane by 6 radiopaque gold markers.

└



ALLEGRA™






Maintaining Coronary Access.

Whatever the anatomic situation, coronary access should be maintained. In all studies conducted with **ALLEGRA™** so far, there is a 0% incidence of coronary obstruction¹⁻³.





In patients with wide coronary sinuses, intra-annular fixation of the ALLEGRA™ and the absence of any interaction with the ascending aorta both ensure continued access to the coronary ostia.

In patients with narrower coronary sinuses, the outflow section of the ALLEGRA™ has large stent cells which further facilitate coronary access. The stent frame can be safely crossed with standard interventional catheters including up to 12 Fr⁶.

In patients with very low-lying coronary ostia, the 6 easily visualized gold markers which mark both the new valve plane and the top of the 12 mm high sealing skirt allow the operator to accurately position the valve taking into account the coronary heights.



Durability by Design.

Durability of a transcatheter heart valve is critically important in younger patients.






The stent frame of the ALLEGRA™ allows pole movement which reduces leaflet stress. Pole movement allows the commissural points of the leaflets to move with every cardiac cycle and reduces stress to the leaflets⁷. In the surgical aortic valve literature reduction of leaflet stress is associated with increased durability⁸.




stent tip deflection during the cardiac cycle




stent frame allows pole movement





For the ALLEGRA™ bovine pericardium is used for all tissue components. Bovine pericardium was chosen as it has an excellent performance profile with very good biocompatibility, outstanding hemodynamic characteristics, low complication rates and extensive data on durability from the surgical literature⁹.

The ALLEGRA™ bovine pericardium is treated to reduce calcification potential. The protective treatment is based on proven techniques such as the elimination of the phospholipid layer of the pericardial cells and the reduction of the glutaraldehyde free bonds. Both of these treatments offer anti-calcification characteristics to the pericardium.



**Bespoke
Delivery System.**

ALLEGRA™ comes with a bespoke delivery system for occlusion free Permaflow implantation technique.

ALLEGRA™






A highly flexible 15 Fr catheter with an 18 Fr valve cartridge allows enhanced cross-ability through the aortic arch.

highly flexible catheter



The ALLEGRA™ delivery system offers an occlusion free Permaflow deployment technique with no hemodynamic compromise during implantation. Permaflow means that at all times during delivery and deployment, flow is maintained providing the clinician with stable hemodynamics and ample time to optimize device position.




occlusion free Permaflow




The ALLEGRA™ delivery system incorporates an easy 3-step deployment. Step 1 is Permaflow, step 2 is the inflow release and step 3 is the final detachment of the valve.



Permaflow




Inflow release



Final detachment of the valve



Visibility for precise implantation is important. The ALLEGRA™ delivery system and the ALLEGRA™ transcatheter heart valve offer unique visibility features which allow precise positioning and implantation of the valve. The delivery system has two ring markers for top and bottom reference and the valve has 6 radiopaque gold markers to indicate the new valve plane and the top of the 12 mm sealing skirt.



(in mm)	ALLEGRA™ 23	ALLEGRA™ 27	ALLEGRA™ 31
Tissue annulus diameter	19 – 22	22 – 25	25 – 28
Annulus perimeter	59.7 – 69.1	69.1 – 78.5	78.5 – 88
Annulus area	280 – 380 mm ²	380 – 490 mm ²	490 – 620 mm ²
Inner diameter aortic valve bioprosthesis	16.5 – 21.5	22 – 24.5	25 – 28

References

1. Wenaweser et al, Transcatheter aortic valve implantation with the NVT Allegra transcatheter heart valve system: first-in-human experience with a novel self-expanding transcatheter heart valve, EuroIntervention 2016
2. Lemos, Clinical performance of a novel transfemoral, supra-annular, early functional, retrievable transcatheter aortic valve system, PCR London Valves 2017 presented at Late Breaking Trials
3. Schaefer et al, Thirty-day outcomes of a novel transcatheter heart valve to treat degenerated surgical valves: the VIVALL multicentre, single-arm, pilot study; EuroIntervention, 2019
4. Data on file.
5. Sathananthan et al, Impact of implant depth on hydrodynamic function with the ALLEGRA transcatheter heart valve following valve-in-valve intervention. EuroIntervention 2019
6. Data on file.
7. Data on file.
8. Christie GW, Barratt-Boyes BG. On stress reduction in bioprosthetic heart valve leaflets by the use of a flexible stent. Journal of cardiac surgery 1991;6:476-81.
9. Yap et al, Aortic valve replacement: is porcine or bovine valve better?, 2013

A collaboration



ALLEGRA™ Delivery System TF

- 15 Fr catheter with an 18 Fr valve cartridge
- Compatible with a 0.035 inch (0.889 mm) guide wire
- Usable length 115 cm
- Hydrophilic coating

Disclaimer

ALLEGRA is a trademark or registered trademark of NVT AG. ALLEGRA™ is CE Mark approved. Other cited trademarks are the property of their respective owners. Data on file at NVT AG for any sustained claims in this document. CAUTION: The law restricts these devices to sale by or on the order of a physician and these products are intended for the use by or under the direction of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Not available in the United States and any other country where applicable health authority product registration has not been obtained. Information contained herein only for presentation outside the US and France. ALLEGRA™ is a product of NVT AG. Blue Sail Medical Co., Ltd is the parent company of NVT AG and Biosensors Intl. and its affiliates are collaborating for the commercialization of the ALLEGRA™ device. © 2021 NVT AG. All rights reserved.

www.nvt-med.com

NVT AG

Rue de Lausanne 29
1110 Morges (Switzerland)
Tel. +41 (0) 21 804 80 41
info@nvt-med.com



NVT GmbH

Lotzenäcker 17
72379 Hechingen (Germany)
Tel. +49 (0) 7471 98979-0
info@nvt-med.com

